

VACUETTE®

Kapillärrör Minicollect



Artikel	Kork/Drag	Volym	Förpackning
450542 Minicollect Z ingen tillsats	Vit	-----	50/500/2000
450534 Minicollect Z serum	Röd	0,5/1,0ml	50/500/2000
450533 Minicollect Z serum gel	Guld	0,5/0,8ml	50/500/2000
450536 Minicollect Lithium Heparin	Grön	0,5ml	50/500/2000
450537 Minicollect Lithium Heparin	Grön	1,0ml	50/500/2000
450535 Minicollect Lithium Heparin gel	Mintgrön	0,8ml	50/500/2000
450530 Minicollect K3 EDTA	Lila	0,25/0,5ml	50/500/2000
450531 Minicollect K3 EDTA	Lila	1,0ml	50/500/2000
450532 Minicollect K2 EDTA	Lila	0,25/0,5ml	50/500/2000
450540 Minicollect FX Natrium Fluorid	Grå	0,25ml	50/500/2000
450541 Minicollect FX Natrium Fluorid	Grå	0,5ml	50/500/2000
450538 Minicollect 9NC Koagulation	Ljusblå	0,5ml	50/500/2000
450539 Minicollect 9NC Koagulation	Ljusblå	1,0ml	50/500/2000
450417 PREMIUM Bärarrör 13x75mm	-----	-----	100/500