

**PRODUCT SPECIFICATIONS****PRODUCT CODE** 901101S**GENERAL DESCRIPTIONS** flocked nasopharyngeal swab 1 pc individually packaged, sterile**DIMENSION**  
STICK DIAMETER: 1.00 mm ± 0.05 mm  
TOTAL LENGTH: 153.00 mm ± 1.00 mm  
FLOCK DIAMETER: 2.50 mm ± 0.05 mm  
FLOCKED LENGTH: 53.00 mm ± 1.00 mm  
BREAK POINT: 100 mm ± 1.00 mm  
STICK MATERIAL: ABS**PACKAGING SPECIFICATIONS****INNER BOX QTY** 300**OUTER BOX QTY** 10**QTY TOTAL** 3000**PACKAGING** PAPER-FILM POUCH**STERILIZATION METHOD** Irradiation**STORAGE TEMPERATURE** 2/30°C - DO NOT FREEZE**SHELF LIFE** 24 MONTHS