

HIGH FLEXIBILITY AND HANDLING SAFETY



HOLDEX®

Single-Use Holder PP

The HOLDEX® Single-Use Holder is a sterile device with a male Luer slip connector to enable compatibility with standard female Luer connections. The excentric position of this Luer connection offers the further advantage of a comfortable, flat puncture angle.

HOLDEX® was especially developed for blood collection from difficult vein conditions. Due to the Luer connector which is positioned off-centre (excentric), the user has a pleasant, flat puncture angle. This increases comfort for both user and patient.

- / An innovative solution for streamlining the product line.
- / One blood collection holder for different puncture devices.
- / Optimal puncture angle for ease of use.
- / Optical blood flow indicator directly at the puncture site.

ONE HOLDER FOR VARYING APPLICATIONS.

- / in combination with safety hypodermic needles
- / in combination with safety blood collection sets

The practical slip connection means easy application with the chosen blood collection safety device.



PRODUCT OVERVIEW

FOR FURTHER INFORMATION AND/OR SAMPLE ORDERING
PLEASE VISIT OUR WEBSITE OR CONTACT US.

HOLDEX® Single-Use Holder PP

Item No.	Description	Inner / Outer [Qty.]
450241	single-packed, sterile, not made with natural rubber latex	100 / 800
450263	single-packed, sterile	100 / 800

RECOMMENDED HANDLING

Please observe the instructions for use at www.gbo.com.



1

The off-centre male connector provides a pleasant and gentle puncture angle.



2

The transparent design allows visual blood flow control directly during puncture.

Note: Please note that the correct venipuncture technique and sufficient venous pressure are necessary conditions for filling the viewing window. The Hypodermic Needle shown here is not sold by Greiner Bio-One.